











Goldent

acrylic teeth mould chart



The 2 layered, cross-linked GOLDENT PREMIUM, and the double cross-linked multi layered GOLDENT EXCLUSIVE acrylic teeth could be ordered in 16 different colour shades according to the standards of vita from type A1-D4, in 44 shapes (21 upper incisors, 9 lower incisors, 7 lower-upper molars)



GOLDENT EXCLUSIVE




Strong surface relief
Light enamel spots
Internal mamelons
Darker necks
Darker canines
Translucent incisal parts
Enamel striation



articulation combination chart

ARTIKULATIONS TABELLE TABLA DE ARTICULACION TABLEAU D'ARTICULATION
 OVERSICHT OVER ARTIKULATION TABELLA DI ARTICOLAZIONE ARTICULATIETABEL
 PARENTATAULUKKO ΠΙΝΑΚΑΣ ΑΡΘΡΩΣΕΩΝ TABEL DE ARTICULAÇÃO KOMBINÁCIÓS TÁBLÁZAT

R1	Z1	Z2	S1
R2	Z2	Z3	S2
R3	Z3	Z4	S1
R4	Z3		S3 M1
R5	Z4		S4 M2
R6	Z4		S3 M1
R7	Z5		S4 M2
R8	Z8	Z9	S5
U1	Z1		S1
U2	Z4	Z3	S2
U3	Z4		S4 M2
U4	Z6	Z5	S4 M2
U5	Z7	Z8	S5
U6	Z8		S5
U7	Z9		S5
X1	Z1		S1
X2	Z6	Z5	S2
X3	Z5		S4 M2
X4	Z7		S4 M2
X5	Z7		S4 M2
X6	Z7		S5



The GOLDENT PREMIUM and GOLDENT EXCLUSIVE acrylic teeth possess CE certification and comply with the requirements of EN ISO 3336 international standard for acrylic teeth, production is performed according to ISO 9001 and ISO 13488 quality standards.